

ClosureFAST™ and ClosureRFS™ Accessory and Supply Ordering Information

Supply Item	Description/Comments	Ordering Information
Procedure Packs		
VNUS Procedure Pack with no introducer sheath set	<ul style="list-style-type: none"> • Sterile surgical supplies • Peel pouch 	Part #: CFP VNUS Medical Technologies 888-797-8346
Introducer Sheaths Sets		
VNUS 6F Introducer Sheath Set, 7cm length	<ul style="list-style-type: none"> • 6F 7cm introducer sheath set w/ .035" 45cm DE guidewire, 19G access needle 	Part #: IS-6F07 VNUS Medical Technologies 888-797-8346
VNUS 7F Introducer Sheath Set, 7cm length	<ul style="list-style-type: none"> • 7F 7cm introducer sheath set w/ .035" 45cm DE guidewire, 19G access needle 	Part #: IS-7F07 VNUS Medical Technologies 888-797-8346
VNUS 8F Introducer Sheath Set, 7cm length	<ul style="list-style-type: none"> • 8F 7cm introducer sheath set w/ .035" 45cm DE guidewire, 19G access needle 	Part #: IS-8F07 VNUS Medical Technologies 888-797-8346
VNUS 6F Introducer Sheath Set, 11cm length	<ul style="list-style-type: none"> • 6F 11cm introducer sheath set w/ .035" 45cm DE guidewire 19G access needle 	Part #: IS-6F11 VNUS Medical Technologies 888-797-8346
VNUS 7F Introducer Sheath Set, 11cm length	<ul style="list-style-type: none"> • 7F 11cm introducer sheath set w/ .035" 45cm DE guidewire, 19G access needle 	Part #: IS-7F11 VNUS Medical Technologies 888-797-8346
VNUS 8F Introducer Sheath Set, 11cm length	<ul style="list-style-type: none"> • 8F 11cm introducer sheath set w/ .035" 45cm DE guidewire, 19G access needle 	Part #: IS-8F11 VNUS Medical Technologies 888-797-8346
Micropuncture Kits		
VNUS 4F Micropuncture Kit	<ul style="list-style-type: none"> • 4F, 7cm micropuncture kit w/ dilator, .018" 45cm guidewire, 21G access needle • Use to size up to a 6F or 8F sheath 	Part #: MPK-4F VNUS Medical Technologies 888-797-8346
VNUS 5F Micropuncture Kit	<ul style="list-style-type: none"> • 5F, 7cm micropuncture kit w/ dilator, .018" 45cm guidewire, 21G access needle • Use to size up to a 6F or 8F sheath 	Part #: MPK-5F VNUS Medical Technologies 888-797-8346

ClosureFAST™ and ClosureRFS™ Accessory and Supply Ordering Information

Supply Item	Description/Comments	Ordering Information
Micro Introducer Set		
VNUS 6F Micro Introducer Set	<ul style="list-style-type: none"> Combines the features of a micropuncture kit and standard introducer sheath 6F, 7cm sheath w/ dilator tapering to 2F, .018" 45cm guidewire, 21G access needle 	Part #: MIS-6F07 VNUS Medical Technologies 888-797-8346
VNUS 7F, 7cm Micro Introducer Set	<ul style="list-style-type: none"> Combines the features of a micropuncture kit and standard introducer sheath 7F, 7cm sheath w/ dilator tapering to 2F, .018" 45cm guidewire, 21G access needle 	Part #: MIS-7F07 VNUS Medical Technologies 888-797-8346
VNUS 7F, 11cm Micro Introducer Set	<ul style="list-style-type: none"> Combines the features of a micropuncture kit and standard introducer sheath 7F, 11cm sheath w/ dilator tapering to 2F, .018" 45cm guidewire, 21G access needle 	Part #: MIS-7F011 VNUS Medical Technologies 888-797-8346
VNUS 8F Micro Introducer Set	<ul style="list-style-type: none"> Combines the features of a micropuncture kit and standard introducer sheath 8F, 7cm sheath w/ dilator tapering to 2F, .018" 45cm guidewire, 21G access needle 	Part #: MIS-8F07 VNUS Medical Technologies 888-797-8346
Guidewires		
VNUS PTFE-Coated Guidewire	<ul style="list-style-type: none"> .025", 150cm length PTFE (Teflon)-coated Double-ended: 1.5mm J and straight Used to help navigate the ClosurePLUS catheter 	Part #: GW-025-150-DE-1.5J VNUS Medical Technologies 888-797-8346
VNUS PTFE-Coated Guidewire	<ul style="list-style-type: none"> .025", 260cm length PTFE (Teflon)-coated Double-ended: 1.5mm J and straight Used to help navigate the ClosurePLUS catheter 	Part#: GW-025-260-DE-1.5J VNUS Medical Technologies 888-797-8346
VNUS PTFE-Coated Guidewire	<ul style="list-style-type: none"> .018, 45cm length PTFE (Teflon)-coated Used to help navigate the ClosurePLEX device 	Part#: GW-018-45 VNUS Medical Technologies 888-797-8346
VNUS PTFE-Coated Guidewire	<ul style="list-style-type: none"> .035, 45cm length PTFE (Teflon)-coated Used to help navigate the ClosureRFS device 	Part#: GW-035-45 VNUS Medical Technologies 888-797-8346
Glidewire		
Boston Scientific	<ul style="list-style-type: none"> 150 and 260cm lengths .025" 	Part #: 630-093 (150 cm) or 5199 (260 cm) 800-225-3226
Y-Connector		
B. Braun Medical, Inc.	<ul style="list-style-type: none"> 9.5F lumen Single Y-connector w/ plastic guidewire insertion tool 	Part #: 610 402 B. Braun Medical, Inc Bethlehem, PA 800-854-6851

ClosureFAST™ and ClosureRFS™ Accessory and Supply Ordering Information

Supply Item	Description/Comments	Ordering Information
Esmark Wraps		
Bioseal	<ul style="list-style-type: none"> • Latex-free • 6" by 4 yards • Sterile 	Part #: 04964/20 Bioseal Placentia, CA 714-528-4695
DeRoyal Industries	<ul style="list-style-type: none"> • Latex-free • 6" by 4 yards • Sterile • CE mark 	Part #: 35-612 20 per case DeRoyal Industries Lake Powell, TN 800-251-9864
Needles		
VNUS Percutaneous Entry Needle	<ul style="list-style-type: none"> • 19G UTW (ultra-thin wall) • Accommodates up to a .035" guidewire 	Part #: AN-19G VNUS Medical Technologies 888-797-8346
VNUS Percutaneous Entry Needle	<ul style="list-style-type: none"> • 19G UTW (ultra-thin wall) with echogenic tip • Accommodates up to a .035" guidewire 	Part #: AN-19GE VNUS Medical Technologies 888-797-8346
VNUS Percutaneous Entry Needle	<ul style="list-style-type: none"> • 21G • Accommodates up to a .018" guidewire 	Part #: AN-21G VNUS Medical Technologies 888-797-8346
VNUS Percutaneous Entry Needle	<ul style="list-style-type: none"> • 21G with echogenic tip • Accommodates up to a .018" guidewire 	Part #: AN-21GE VNUS Medical Technologies 888-797-8346
12G IV catheter	<ul style="list-style-type: none"> • 12G Angiocath - Becton Dickenson • For use with ClosureRFS vein access 	Progressive Medical Inc 800-764-0636 order # 18556 Blue Ridge Medical 800-554-9587 order # BD38-2830-2 Moore Medical Supplies 800-234-1464 order # 71942
Transducer Cover		
VNUS Probe Cover	<ul style="list-style-type: none"> • 6" x 96" • Contains sterile gel packet and rubber bands • Latex-free 	Part #: PC-6x96 VNUS Medical Technologies 888-797-8346

ClosureFAST™ and ClosureRFS™ Accessory and Supply Ordering Information

Supply Item	Description/Comments	Ordering Information
Ultrasound Gel		
Sterile Gel Packets	<ul style="list-style-type: none"> • 20 ml packets, 48/box 	Part #: 610-117 Civco Medical Instruments Kalona, Iowa 800-445-6741
Non-sterile gel	<ul style="list-style-type: none"> • .25 liter bottles, 12/case • 5 liters w/ 1.25 liter dispenser bottle 	Part #: 610-130 (.25L bottles) 610-116 (5L dispenser bottle) Civco Kalona, IA 800-445-6741
Pressure Bag		
1000ml IV Pressure Infusor Bag	<ul style="list-style-type: none"> • With pressure gauge 	Part #: 7789 Alliance Medical Russellville, MO 888-633-6908
500ml IV Pressure Infusor Bag	<ul style="list-style-type: none"> • With pressure gauge 	Part #: 7788 Alliance Medical Russellville, MO 888-633-6908